


## Lab 11: Electrochemical Cells

Name \_\_\_\_\_

Oxidizing Partner's name \_\_\_\_\_

**Purpose:** To build an electrochemical cell with Zn, Cu and the appropriate solutions, to measure its voltage and to find ways of increasing it.

### Procedure:



1. Add 0.5 M  $\text{CuSO}_{4(aq)}$  or  $\text{Cu}(\text{NO}_3)_{2(aq)}$  solution to a 100 mL beaker until it is half-full. Place a clean strip of copper into the solution, and connect it to the positive end of the voltmeter.
2. Take a second beaker and fill half of it with  $\text{Zn}(\text{NO}_3)_2$  solution. Place a clean strip of zinc into the solution, and connect it to the negative end of the voltmeter.
3. Get the salt bridge ready by filling a U-tube with an electrolyte. Seal both ends with cotton, and try to minimize the amount of air in the tube.
4. Record the voltage.


Voltage reading(use one estimated figure; include unit)	
---	--

5. On the next page, record any changes observed in the next 5 minutes. Try lifting the salt bridge out of the solution, and then placing it back in.


**Lab 11: Electrochemical Cells**

Name \_\_\_\_\_

Oxidizing Partner's name \_\_\_\_\_



6. Use the multi meter to record current \_\_\_\_\_
7. If you noticed no changes yet, why do you think that is? \_\_\_\_\_
8. Try connecting your cell and that of your neighbour's *in a way that increases voltage*. What voltage is recorded? \_\_\_\_\_
9. Draw the setup of your battery (the way it's connected to your neighbour's.) Technically that's what a battery is: two or more connected cells.



**Lab 11: Electrochemical Cells**

Name \_\_\_\_\_

Oxidizing Partner's name \_\_\_\_\_

**Analysis:**

1. a. What metal served as the anode? \_\_\_\_\_  
b. What in the experiment revealed what the anode was? \_\_\_\_\_

2. Write the two half-reactions that occurred:

\_\_\_\_\_


\_\_\_\_\_

3. Write the overall balanced reaction.

\_\_\_\_\_

4. Was there any voltage before inserting the salt bridge? \_\_\_\_\_

5. Use the measurement of current flowing through your electrochemical cell. Keeping in mind that exactly  $1A = 1C/s$ , and that **every mole of electrons** has a charge of 96 400 C, how many hours would you have to wait before you saw 1.0 g of copper at the cathode?



**An Ox=**

**anode = oxidation**

**Conclusion:**